



Introducing the

# KaVo SURGtorque™

A new angle that's distinctively KaVo



The new KaVo SURGtorque™ features 19W of cutting power for faster preps, side and rear exhaust to help increase patient safety, and enduring KaVo quality.

With its uniquely angled head, increased access is simple. Experience for yourself, the difference between ordinary and extraordinary.



Product	Description	Part No.
SURGtorque S 459 L	45 degree head angle, 19 W, 2 side vents, ceramic bearings, single-port spray, fiber optic MULTIflex LUX connection	3.000.4808
SURGtorque S 459 C	45 degree head angle, 19 W, 2 side vents, ceramic bearings, single-port spray, 4-hole Midwest connection	3.000.5063




Precise Power  
 Uniquely Angled Head  
 Increased Access



KaVo. Dental Excellence.



# A New Angle That's Distinctively KaVo



## POWER & PRECISION

### PRECISE POWER

#### Reliable performance and durability

- 19 W @ 2.8bar
- Water Jet (not Mist) for perfect bur cooling
- Ceramic bearings – for long lifetime and lower sound level
- Precise work with a balanced instrument, carbide lined chuck for more concentric cutting and more precise preps
- MADE IN GERMANY
- 12 Month warranty
- 5 surface coating layers – prepared for over 1000 sterilization cycles

## 45° HEAD

### UNIQUELY ANGLED HEAD

#### Ergonomics Handling and ease of use

- 45° Head design for perfect access to the molars compatible for MULTIflex / 4-hole standard
- Light version for perfect visibility
- Low weight, ergonomic balanced for longer treatments

## SAFETY & HYGIENE

- Exhaust Air diverted from cutting area – No air contamination into cutting area
- Smooth surfaces free of recess – easy to clean and hygienic design
- 5 surface coating layers – prepared for over 1000 sterilization cycles
- Up to 30N chuck holding force – strong and safe chuck system, perfect for a longer bur