




Professional Digital Dentistry  
DLP 3D Printer for Dental

# D2




Distinct Advantage

“ Qualified Accuracy and Repeatability ”



D2 is manufactured using digital calibration systems, all of which are individually manufactured to ensure the highest level of accuracy.

Outstanding Features




## Material & Application


### | Working Model for Digital Impression (DM-1)<sup>gray</sup>



### | Temporary (DT-1)



### | Cast Pattern (DC-1)




### | Set-up Model for Clear Aligner (DM-1)<sup>yellow</sup>



### | Denture Base (DD-1)



### | Splint for Anti-snoring guide (DS-1)



### | Implant Guide (DG-1)



※ All these dental prostheses were printed on a D2 using materials developed by Veltz3D.


## D2 Specifications

Specification	D2-120	D2-150
Layer thickness	25 ~ 100 $\mu$ m	
Build size	120(W) X 67.5(D) X 100(H) [mm]	150(W) X 84.3(D) X 100(H) [mm]
Printing speed	Up to 25mm/h, Z-axis @ 50 $\mu$ m layer thickness	
Optical engine	1920 X 1080 (Full HD), 405nm UV LED	
XY resolution	62.5 $\mu$ m	78 $\mu$ m
Power	100~240 Vac, 50/60Hz / Approx. 150W	
Dimension/Wt.	310(W) X 430(D) X 910(H) mm / Approx. 40kg	




## MP100 / Post Curing Machine (UV curing light)

Applicable resin for FDA Class II (Perfect curing performance with wide enough wave length)  
 Guarantee the quality through optimized mode selection of each specific resin material



▲ soft touch switch



▲ mode information

Item	Post Curing Machine (MP100)	Light	Metal UV lamp
Size	220 x 380 x 150[mm]	Tray size	173 x 140 x 25[mm]
Power	220vac. 50-60Hz.	Weight	6.15kg
Current	0.34 A	Watt	74 W
Acc.	Tray, User's manual, Power cable		



# HEPHZIBAH

## CONTACT US

sales@veltz3d.com  
 (22121)86, Gilpa-ro 71beon-gil, Nam-gu, Incheon, Korea  
 Tel. +82-32-509-5822 Fax. +82-32-502-5519

Copyright © 2018 Hephzibah co., Ltd. All Rights Reserved.

<https://stomshop.pro>

